NeuRx Diaphragm Pacing System™
Clinician's Guide
Système de stimulation NeuRx™
Guide du clinicien
NeuRx™ System zur Atemstimulation
Handbuch für den Arzt
Sistema di stimolazione NeuRx™
Guida per il medico
Sistema de estimulación NeuRx™
Guía para médicos
Stimulační systém bránice NeuRx™
Příručka pro lékaře

TABLE OF CONTENTS

1.0	INTENDED USE.....	3
2.0	EXPLANATION OF SYMBOLS	3
3.0	WARNINGS.....	6
4.0	CAUTIONS.....	8
5.0	DEVICE DESCRIPTION AND PRINCIPLES OF OPERATION	8
6.0	TECHNICAL AND PHYSICAL DESCRIPTION.....	10
7.0	IMPLANTABLE COMPONENTS	10
8.0	EXTERNAL COMPONENTS.....	12
9.0	CLINICAL STATION.....	14
10.0	SURGICAL COMPONENTS	15
11.0	INFORMATION	18
12.0	PROCEDURE RISKS.....	18
13.0	TRAINING OF SURGICAL IMPLANTATION.....	19
14.0	DESCRIPTION OF LAPAROSCOPIC SURGICAL PROCEDURE	19
15.0	MAPPING OF THE DIAPHRAGM	20
16.0	INSERTION AND PLACEMENT OF THE ELECTRODES	21
17.0	ROUTING.....	25
18.0	BLOCKING	25
19.0	APPLY STRAIN RELIEF BOOT	31
20.0	REFERENCES	35

Additional instructions delivered with the NeuRx DPSⁱ that physicians should be familiar with:


PN 77-0035-EN Patient / Caregiver Instruction Manual

1.0 INTENDED USE


The NeuRx Diaphragm Pacing System (NeuRx DPS[®]) provides ventilatory support in patients with diaphragm dysfunction of neuromuscular origin who require chronic ventilatory support because of diaphragm paralysis. The device is indicated in patients with high-level spinal cord injury, whose remaining phrenic nerve, lung and diaphragm muscle function are sufficient to accommodate electrical stimulation.

2.0 EXPLANATION OF SYMBOLS

-  The **Warning** symbol precedes warning information that mitigates a risk that is not obvious to the operator. Indicates that a potentially hazardous situation which, if not avoided, could result in harm to the operator or patient.
-  The **Caution** symbol appears next to precautionary information when the intention is solely to inform. Indicates that a potentially hazardous situation which, if not avoided, may result in minor or moderate personal injury or property damage. This word is used to also alert against unsafe practices.
-  The **Manufacturer** symbol appears next to the manufacturer's name and address.
-  The **Reference** symbol appears preceding the part number for the device. The part number is a unique numeric identifier for the device.
-  The **Lot** symbol appears preceding the lot number for a device. Devices manufactured at the same time using identical material and parts will share a common lot number.
-  The **Serial Number** symbol appears on devices that require unique identification.
-  The **Use Until** symbol appears on devices that have an indication of the date by which the device should be used. The date is expressed as the year and month, with the month referring to the end of the month.
-  The **Manufactured Date** symbol appears on devices as an indication of the date of manufacture. The date is expressed as the year and month.
-  The **Temperature Limits** symbol appears on packages of devices as an indication of the storage and transit temperature limits. While the symbol appears on actual devices as an indication of the operational temperature limits.

-  The **Keep Dry** symbol appears on all packages of devices requiring to protect the packaging from potential damage.
-  The **Don't Use If Packing Damaged** symbol appears on all packages of devices requiring to dispose of the device if the packaging has suffered damage.
-  The **Accompanying Documents** symbol appears on all packages of devices indicating that instructions for use are available for additional information.
-  **CE₂₇₉₇** The **Regulatory Marking of Conformity** symbol indicates that the device meets Medical Device Directive 93/42/EEC. This has been certified by notified body number 2797.
-  The **European Community Representative** symbol indicates the identification of the authorized representative for the distribution of devices into the European community.
-  The **Type BF Applied Part** symbol appears on powered equipment that connects directly to a patient. It is an indication of the degree of protection provided against electric shock, patient leakage current and patient auxiliary current.
-  The **Warning** symbol on powered equipment indicates physiological effects not obvious to the user that can cause harm.
-  The **On / Off** symbol on powered equipment indicates push-button ON/OFF power control of the device.
-  The **Consult Accompanying Documents** symbol appears on powered equipment indicating that instructions for use must be consulted for safety.
-  **IP₂₄** The **Ingress Protection (IP) Classification** symbol appears on powered equipment indicating that the device is protected from splashing water.
-  MR Unsafe. A device that is known to pose hazards in all MR environments.
-  Gamma Sterilization. This simple, proven process is safe, reliable, and highly effective at treating single-use medical devices. With the ability to penetrate

products while sealed in their final packaging, gamma irradiation supports the manufacturing and distribution process by facilitating final packaged product as well as raw materials, whilst still ensuring full sterility of the product.

Indicates a medical device that has been sterilized using ethylene oxide.

3.0 WARNINGS

- This device should be used only under the direction of a physician.
- The patient must always have a back-up means of breathing help (ventilation) available. If the patient thinks there is any problem with the NeuRx DPS, then they should use their back-up means right away. They should turn off the NeuRx DPS when they switch to their back-up means of breathing help.
- NeuRx DPS could interfere with some medical equipment. Some medical equipment could interfere with the NeuRx DPS:
 - **All active implantable medical devices.** This will include devices such as implanted cardiac pacemakers, implanted cardioverter defibrillators (ICDs), implanted neurostimulators, and body worn medical devices (e.g., insulin pump). Use of the NeuRx EPG stimulator may interfere with these devices.
 - **Surgery.** Use of high-frequency surgical equipment may cause burns where the electrode wires pass through the skin. It might also damage the NeuRx EPG if connected.
 - **Diathermy treatment.** Diathermy treatment should not be performed within 30cm of the implanted electrode leads. Unwanted tissue heating through the electrode leads could occur.
 - **External electrical stimulation** such as **transcutaneous electrical nerve stimulation (TENS)**. Such stimulation should not be done in the chest area near the electrode wires. Unwanted diaphragm contraction could occur.
 - **Shortwave or microwave therapy.** Operating the NeuRx DPS close to (about 3 feet from) such equipment may interfere with your NeuRx DPS.
 - **MR Magnetic Resonance Imaging (MRI) test.** The PermaLoc® electrode is MR Unsafe. Do not perform a MRI test while implanted with the PermaLoc electrodes.
 - **MR Magnetic Resonance Imaging (MRI) test.** The NeuRx EPG stimulator and surface electrodes are MR Unsafe. The NeuRx DPS has not been tested with MRI. MRI could cause the electrode wires to move. MRI could also cause unwanted tissue heating through the electrode wires.
- Avoid trans-thoracic stimulation.
- Avoid accidental contact between connected but unused applied parts (cable or leads) and other conductive parts including those connected to protective earth.
- The long-term effects of chronic electrical stimulation are unknown.

WARNINGS (contd)

- Safety has not been established for the use of the device during pregnancy.

- It should not be used in patients with suspected heart problems or epilepsy.
- It should not be used if there are skin eruptions such as phlebitis, thrombo phlebitis, or varicose veins.
- **ELECTROMAGNETIC INTERFERENCE WARNING:** Some electrically powered equipment gives off electromagnetic waves which could interfere with your NeuRx EPG. When using your NeuRx EPG around electrical equipment, check the NeuRx EPG screen to make sure the EPG is working.
- Avoid electro-cauterization in the area of the implanted electrodes.
- Disconnect the NeuRx EPG during electrical diagnostic tests such as EMG or ECG.
- Do follow the electromagnetic compatibility (EMC) information provided. The NeuRx EPG needs special precautions regarding EMC. To reduce the possibility of interference on the NeuRx EPG from other electrical equipment or the NeuRx EPG affecting other electrical equipment, do NOT use cables or accessories with your NeuRx EPG other than those specified.
- Do not use this device if skin in the area is swollen, infected, or inflamed.
- **RF COMMUNICATION WARNING:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the NeuRx EPG stimulator, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- **FLAMMABILITY WARNING:** Do NOT use the NeuRx EPG in an oxygen enriched environment, such as a hyperbaric oxygen chamber or near a flammable anesthetic mixture with air, oxygen or nitrous oxide. The NeuRx EPG is not categorized as AP (anesthetic-proof) or APG (anesthetic-proof category G - gas) type of equipment.
- Do NOT operate the NeuRx EPG stimulator within 30cm of security systems, metal detectors, or Electronic Article Surveillance (EAS).
- The patient should avoid eating or drinking while conditioning with the NeuRx DPS. There is a risk of food or liquid entering their lungs. If the patient uses the NeuRx DPS full time, discuss ways to reduce the risk of aspiration.
- Visually inspect all STERILE barriers before use. Do NOT use the device if the sterile barrier is open, damaged or broken.

3.1 CONTRAINDICATIONS

The following are conditions contraindicating NeuRx DPS use in patients that:

The device is not intended to be used in patients who do not have a viable phrenic nerve

The device is not intended to be used in pregnant women

4.0 CAUTIONS


- This device should be kept out of the reach of children, pets or pests.
- The NeuRx DPS is a body worn device kept close to the body and the patient cable should NOT be allowed to hang loosely. This is to avoid the potential of any strangulation. The patient cable should be positioned appropriately to not allow any hanging loops.
- Do not expose the pacing device to excessive moisture, heat, extended sun light, excessive dust or lint, or severe mechanical shock. If display indicates system failure, pain is felt at the electrode exit site, or device exposed to excessive moisture, heat or shock, the patient should use the back-up means of breathing help and stop use of the device.
- Precautions should be observed when there is a tendency to hemorrhage following acute trauma or fracture, following recent surgical procedures when muscle contraction may disrupt the healing process, or where sensory nerve damage is present.
- Some patients may feel skin irritation, sensitivity, or have an allergic reaction because of the NeuRx DPS.

This may be due to:

- The stimulation
- The adhesive on the skin bandage.
- The transparent dressing used over the gauze that covers the electrode wires. (Tegaderm[®] and Op-Site[®] are examples of transparent dressings). Irritation can usually be reduced by changing the settings of the EPG or removing the adhesive.
- To protect the NeuRx EPG against damage due to moisture, **Do NOT allow the NeuRx EPG to get wet.** This includes during bathing, showering, swimming, or any other activity in which you could get wet. The EPG may quit working and not be available for use when needed.
- If the pacing device is not to be used for an extended period or in storage, please remove the primary battery to prevent any possible damage.

5.0 DEVICE DESCRIPTION AND PRINCIPLES OF OPERATION

- The NeuRx Diaphragm Pacing System[®] is a system designed to help patients breathe by using electrical stimulation to contract the muscles.
- It is implanted using standard laparoscopic surgical techniques in an outpatient procedure
- The implanted intramuscular diaphragm electrodes are connected to the NeuRx[®] External Pulse Generator (EPG) (external) using the patient cable and the electrode connector, which is secured to the patient using the connector holder bandage.

- The stimulator provides repetitive electrical stimulation to the implanted electrodes to cause the patient's diaphragm to contract and cause the patient to inhale in a manner similar to natural breathing.
- The EPG will be programmed so that it produces the right stimulation patterns. If the stimulation is uncomfortable, adjust the EPG settings to reduce or eliminate the discomfort.
- The user simply connects the device to the implanted electrodes and turns it on for use; no other controls are available or necessary for operation.
- During use, the stimulator should be kept close to the patient's body to avoid pulling on the cable and electrodes. The stimulator can be placed in a pocket of the patient's clothing, a fanny pack or simply placed on a table or other convenient location.

6.0 TECHNICAL AND PHYSICAL DESCRIPTION

The NeuRx Diaphragm Pacing System (NeuRx DPS) is an intramuscular, percutaneous, motor point diaphragm stimulation system. It is implanted using standard laparoscopic surgical techniques. The implanted intramuscular diaphragm electrodes are connected to a four-channel external pulse generator at a percutaneous exit site. The pulse generator provides a capacitively coupled, charge-balanced, biphasic stimulation to each electrode with a common indifferent electrode that is placed subcutaneously. The pulse generator controls the charge delivered through clinician programmed parameters of pulse amplitude, pulse duration, pulse frequency, pulse ramp, inspiration time, and respiratory rate. The clinician uses a clinical station to evaluate electrode response to stimulation and program the external pulse generator with the patient-specific parameters. The user simply connects the device and turns it on for use; no other controls are available or necessary for operation.

The locations for implantation are identified by electrically mapping the inferior aspect of the diaphragm. Using the laparoscopic electrode delivery instrument, the PermaLoc® intramuscular electrodes are surgically implanted in the diaphragm muscle in proximity to branches of the phrenic nerve without making contact or manipulating the nerve. The PermaLoc electrodes are tunneled, subcutaneously, to a percutaneous exit site on the lateral chest region. An indifferent return electrode (anode) is placed subcutaneously and exits at the same chest location. These electrodes are connected to an external pulse generator (EPG) (stimulator) that controls the timing and level of diaphragm pacing stimulation. Each electrode may be controlled individually in terms of charge (pulse duration and pulse amplitude) delivered and grouped together to recruit the diaphragm muscle to elicit the desired level of inspiratory effort.


7.0 IMPLANTABLE COMPONENTS

7.1 Intramuscular Electrode (PermaLoc® electrode)

The PermaLoc intramuscular electrode is a double helix wound lead with exposed 316LVM stainless steel stimulating surface and polypropylene reinforced core. The PermaLoc also has a barb at the implanted end composed of 14 pieces of polypropylene suture fused together. The body of the lead is insulated with PFA (perfluoroalkoxy) fluoropolymer coating and terminated in a 316L stainless steel pin with a silicon reinforcing sleeve. All of the materials have a history of long-term implantable use as part of previously approved products. The lead is 59 - 63cm in overall length, 0.75mm in diameter, and has a 9mm de-insulated stimulating tip. The surface area of the stimulating tip is a minimum of 18.5mm². The electrodes are packaged in Tyvek packaging and sterilized by ethylene oxide exposure.

7.2 Indifferent Electrode (Anode)

The indifferent electrode provides a common return current path for all of the electrodes implanted in the diaphragm. It is implanted in the subcutaneous tissue of the lateral chest region and is tunneled to the percutaneous exit site. The lead is fabricated of the same double helix wound 316LVM stainless steel as the intramuscular electrode and percutaneous extension lead. It is also terminated in a 316L stainless steel pin at one end and 7cm de-insulated for the return electrode tip. The overall length is 19cm with a 0.75mm diameter. The indifferent electrode is packaged in Tyvek packaging and sterilized by ethylene oxide exposure.

8.0 EXTERNAL COMPONENTS

8.1 NeuRx DPS Stimulator

The patient-external pulse generator (EPG) is an external four-channel battery-powered device that controls the stimulus output and respiratory timing. The four output channels are independently controlled, capacitively-coupled, biphasic outputs with a common return. The device is packaged in an impact-resistant plastic enclosure with patient cable connector on the top, display and power buttons on the front and replaceable battery compartment on the back. A programming connector is located in the battery compartment for connection to the clinical station. The EPG has no controls that allow modification to any parameter settings.

8.2 Patient Cable

A five-conductor cable is provided that connects from the external pulse generator to the electrode connector socket. The pulse generator end of the cable is a positive locking ODU medical-grade push-pull connector. The cable is a silicone jacketed, multi-conductor, shielded cable that is 0.75m in length. The electrode connector end is a custom-molded five conductor connector. It is a precision mated strip of pins embedded in the molded PVC. The patient cable conducts the electrical signals from the controller unit to the implanted electrodes. The cable is provided non-sterile to the end-user and the cable assembly has passed biocompatibility testing.

8.3 Electrode Connector Socket Kit

Each electrode is terminated into a precision ITT Canon socket. The individual sockets are crimped onto the electrode and inserted into a carrier strip. The sockets are inserted into the carrier strip in a set sequence to mate with the patient cable. The crimp connections are protected with a strain-relief boot.

8.4 Strain Relief Boot



The boot is a urethane strain relief for the electrodes that connect at the connector block. It is composed of a medical-grade urethane and is molded by Synapse. Once the electrodes are terminated in the connector socket block, the boot shell is backfilled with silicone adhesive. Once the adhesive is filled in the boot, the boot/assembly block is placed into the connector holder. The urethane boot is packaged as part of the Surgical Connector Kit in a metalized barrier bag with a desiccant so as to minimize exposure to moisture. It is supplied non-sterile to the end-user.

8.5 Connector Holder

The connector holder is an adhesive bandage that has a clip on the surface that secures the connection between the electrodes and the patient cable. The adhesive bandage is MED5322 hypoallergenic fabric medical tape intended for applications to the skin for sustained periods. The bandage can remain in place on the skin for up to 7 days.

8.6 Backup Indifferent Electrode Interconnect

A interconnect assembly is provided that allows the temporary connection of a surface indifferent electrode. With redundancy in all other components, this assembly provides redundancy for the implanted anode. In case of a malfunction with the implanted anode, breakage or extraction, this component can be used to replace that function and thereby allow continued therapy without the need for immediate re-insertion of an implanted anode. The availability of a backup indifferent electrode minimizes the interruption of stimulation.

 <p>NeuRx DPS[†] External Pulse Generator (EPG)</p>	 <p>Patient Cables</p>
 <p>Electrode Connector Block</p> <p>The diagram shows a rectangular connector block with four pairs of pins labeled L, L, R, R, and A. Above the pins, there are four small electrodes. Arrows point from each electrode to its corresponding pin. The labels L, L, R, R, and A are positioned below the pins.</p>	 <p>Molded Strain-Relief Boot and Assembly Block Placed in Connector Holder</p>
 <p>Connector Holder</p>	 <p>Backup Indifferent Electrode Interconnect</p>

9.0 CLINICAL STATION

There are three primary aspects of the device implementation that the clinical station provides. The station provides intra-operative mapping functionality and electrode characterization, incorporates DPS stimulator functionality, and external DPS stimulator programming capability.

The first application of the station in the device implementation is to provide intra-operative stimulation and sensing of stimulated response. This surgical mapping mode utilizes the surgical components listed below to provide twitch or burst stimulation to record and display the abdominal pressure response through the solid state pressure sensor. A stimulator mode is used to test the channels individually and in combination at the end of the surgery to make sure that all electrodes are intact and providing the anticipated response.

The station is then used to characterize the response from each electrode utilizing the four stimulus output channels. The amplitude, pulse width, frequency, and pulse ramping response are each characterized to optimize the tidal volume and patient comfort on a per breath basis. The inspiration time and respiratory rate are set for appropriate minute ventilation. The stimulator mode can be used to fine tune the settings for the patient.

Finally, the station sends the desired stimulus and respiratory timing parameters to the external stimulator in the programming mode. The station can also read the currently programmed parameters from an external stimulator to verify the settings.

Reference the Synapse Biomedical Clinician Station Manual (PN 77-0048) for instructions on operation of the Clinical Station.

Clinical Station

10. SURGICAL COMPONENTS

10.1 Sterilization

With the exception of the Clinical Station and Electrode Delivery Instrument, surgical components are provided in a double-pouch packaged in medical-grade Dupont[®] Tyvek[®] pouches and EtO sterilized. The Electrode Delivery Instrument is steam-sterilized at the surgical site and the Clinical Station is not sterilized. The Disposable Electrode Delivery Tool is Gamma (Sterile R) sterilized and packaged in a sealed clear plastic Tyvek sterile tray.

10.2 Mapping Instrument

The initial step in the surgical implementation is the laparoscopic mapping of the diaphragm. This may be performed by introducing and connecting to an available laparoscopic dissector for stimulation or using the optional 5mm mapping instrument. Either instrument is used to stimulate the inferior surface of the diaphragm in a grid pattern to identify optimal implantation sites of the intramuscular electrodes. The connected laparoscopic dissector may be applied to sequential sites on the diaphragm by the surgeon and stimulated. Optionally, the mapping instrument may be applied to sequential sites on the diaphragm by the surgeon and secured by applying the operating room vacuum through the central lumen of the probe. Stimulation is applied in either a twitch or burst mode from the clinical station to elicit an abdominal pressure change.

10.3 Transducer to Trocar Pressure Tube

A one-meter section of PVC tubing, with male Luer lock connectors on either end, is used to connect a Trocar port to the solid-state pressure sensor. The long length of the tube permits connection to the pressure sensor outside of the sterile field. The transducer to trocar pressure tube is packaged in Tyvek packaging and sterilized by ethylene oxide exposure.

10.4 Solid-State Pressure Sensor

A differential, 1 PSI full scale, pressure sensor transduces the abdominal pressure changes to an electrical signal for the clinical station. It connects to the pressure tube with a female Luer lock and to the clinical station with a positive locking medical-grade connector. The electrical signal provides an indication of relative pressure change.

10.5 Cable Set

A set of cables with touch-proof connectors are used to connect off the sterile field from the mapping instrument to the clinical station. A set of 3m meter cables connect to the mapping instrument or clip leads to test implanted electrodes. Another cable connects from the surface anode to the clinical station.

10.6 Surface Anode

The surface anode is an adhesive electrode that is placed on the skin during the intraoperative procedure. The anode is manufactured by Axelgaard and uses a proprietary hydrogel adhesive to adhere to the surface anode to the skin.

10.7 Electrode Delivery Instrument

The Electrode Delivery Instrument is used intraoperatively to implant the electrodes.

Note: **Section 10.7.1** describes the Reusable (after cleaning and re-sterilization by the Clinical site) Electrode Delivery Instrument (23-0031 Ref. IFU 77-0052).

Section 10.7.2 describes the Disposable Electrode Delivery Tool (21-0039-99 Ref. IFU 77-0086). Both tools have the same indications for use.

10.7.1 Reusable Electrode Delivery Instrument

The Reusable Electrode Delivery Instrument is an 11mm durable laparoscopic instrument device. It is manufactured from surgical stainless steel (304SS) which has been coated with chromium using the MEDCOAT 2000 electroplating process. The proprietary MEDCOAT 2000 process applies chromium to the surface of the stainless steel parts creating a low friction surface that is lubricious and biocompatible.

The barbed intramuscular electrode is loaded in the lumen of the instrument with the de-insulated barb extending out of the needle. The skirt of the polypropylene barb is loaded inside of the needle. When the needle is extended and inserted between the muscle fibers, parallel to the diaphragm surface, the de-insulated barb catches on the fibers and the lead is drawn out of the lumen as the instrument is withdrawn. The instrument is provided non-sterile to the end-user with instructions for cleaning and sterilization.

Figure 1: **Reusable Electrode Delivery Instrument**

10.7.2 Disposable Electrode Delivery Tool

The Reusable Electrode Delivery Instrument is a 10mm single patient use disposable laparoscopic instrument. The barbed intramuscular electrode is loaded in the lumen of the instrument with the de-insulated barb

extending out of the needle. The skirt of the polypropylene barb is loaded inside of the needle. When the needle is extended and inserted between the muscle fibers, parallel to the diaphragm surface, the de-insulated barb catches on the fibers and the lead is drawn out of the lumen as the instrument is withdrawn.

Figure 2: Disposable Electrode Delivery Instrument

10.8 Tunnelers

The lead tunnelers are used intraoperatively to guide the electrode to the implantation site. They are thin-walled tubes of stainless steel (304SS) which conform the ASTM A908 standard for needle tubing. As the material composition of the lead tunnelers is similar to the mapping probe cannula, the testing conducted on the mapping probe by Synapse, namely cytotoxicity, sensitization, irritation, and systemic toxicity (acute systemic and pyrogenicity) is considered relevant and supports the biocompatibility of the lead tunnelers. The lead tunnelers are packaged in Tyvek packaging and sterilized by ethylene oxide exposure.

10.9 Clinical Station

The clinical station provides intra-operative mapping functionality, incorporates NeuRx DPS External Pulse Generator functionality, and NeuRx DPS External Pulse Generator programming capability.

10.10 Clinician Crimp Tool

The electrode leads are terminated prior to installation into the connector socket block by using the Clinician Crimp Tool. The termination is of the electrode lead requires the crimping of a larger contact socket that is compatible with the connector socket block.

10.11 Socket Pusher

After the terminated electrode leads after terminated, the Socket Pusher is used to insert each terminated lead into the connector socket block.

10.12 Device Packaging

The components of the NeuRx DPS are kitted into configurations that satisfy end-user requirements. The NeuRx DPS System Kit as delivered contains a Sterile Surgical Kit, Patient Kit and a Destination Kit. The Sterile Surgical Kit includes the implanted electrodes and termination connector parts in addition to all surgical disposable items that are used during the surgical installation. The Patient Kit is given to the patient and includes the external stimulators and the consumable components such as cables, batteries, and connector holders. The Destination Kit contains Patient & Caregiver Instruction and other Patient Information along with a patient survey and warranty information.

The surgical components are delivered sterilized and delivered in DuPont Tyvek® pouches. The Disposable Electrode Delivery Tool is Gamma (Sterile R) sterilized and packaged in a sealed clear plastic Tyvek sterile tray. Visually inspect all STERILE barriers before use. Do NOT use the device if the sterile barrier is open, damaged or broken.

11.0 INFORMATION

See the document %Synapse Biomedical NeuRx Diaphragm Pacing System® Patient Information+(PN 77-0038) provided with the device, which is meant to be provided to the patient.

12.0 PROCEDURE RISKS

- There is a risk of diaphragm penetration during the procedure, which could cause a condition known as capnothorax
- There is a risk of infection and/or inflamed tissue at the electrode implantation sites
- There is a risk of bleeding at the electrode implantation sites
- There is a risk of nerve, tissue or organ damage as a result of the procedure
- There is a risk that the electrode wires could break off in the body leading to reduced or intermittent diaphragm pacing or failure of the pacing system
- There is a risk of cardiac arrhythmia being caused by the placement of the electrodes in the chest cavity
- There is a risk of skin irritation or hypersensitivity from the electrical stimulation or from the tape used with the electrodes or from the skin bandage that holds the electrode connections
- There is a risk that the body may not be compatible with the materials used in the electrodes and their wires
- There is a risk of choking during eating and of sleep apnea if a Passy-Muir® valve is not used during the training period
- There is a risk of discomfort during the conditioning period.
- There is a risk that a patient may not have sufficient muscle reaction when using the stimulator and the product may not work for every patient

- At this time, there is insufficient clinical data to determine safety in implanting patients with cardiac pacemakers, therefore, patients should not be implanted with this device if they have a cardiac pacemaker or other implanted electrical device
- This product should not be used by patients with suspected or real heart problems or who have epilepsy
- The safety of this device in use during pregnancy is unknown
- The long-term effects of electrical stimulation of the diaphragm are unknown
- It is possible that stimulation from the diaphragm pacing system could stop either due to electrode breakage, cable disconnection, or stimulator failure. If one of these happens, breathing will stop. Without prompt attention, this could result in permanent disability or death. This risk is reduced by using back-up electrodes and sounding an alarm whenever the stimulator detects improper operating conditions.
- There is a risk of aspiration when using the device. While becoming used to the stimulation and timing, it is recommended that a one-way speaking valve on the tracheostomy be used while eating or drinking when on the pacing device. It is also recommended to use the one-way speaking valve while sleeping to avoid upper airway obstructions.
- Patients may experience increased spasms with the stimulation while their body becomes used to the stimulation. This typically subsides within a few days of use.

13.0 TRAINING OF SURGICAL IMPLANTATION

New surgeons will be proctored in the implantation of the NeuRx DPS® by a qualified, trained surgeon. To be considered for training, the new surgeon must be credentialed for laparoscopic surgery at their institution. The training procedure will, minimally consist of the qualified trainer proctoring from outside of the surgical field, in the trainee's operating room. The trainer will review the procedure, using this manual, and the patient information to assure that appropriate candidate selection has been made. Specific precautions will be reviewed, including the identification of a capnothorax with preventive measures and steps to resolve, prior to the procedure. Any peri-operative morbidity, mortality, or device malfunctions will be reported according to the FDA Medical Device Reporting (MDR) procedures and reviewed by Synapse Biomedical. Upon completion of training the trainee may request additional onsite proctoring or the trainer may identify that additional proctoring is necessary. The training will be documented on a Surgeon Training Memo, to be kept on record at Synapse Biomedical.

14.0 DESCRIPTION OF LAPAROSCOPIC SURGICAL PROCEDURE

WARNING: Visually inspect all STERILE barriers before use. Do NOT use the device if the sterile barrier is open, damaged or broken.

14.1 INCLUSION PRE-SCREENING

Candidates must have diaphragm dysfunction of neuromuscular origin who require chronic ventilatory support because of diaphragm paralysis. The device is indicated in patients, such as high-level spinal cord injury, whose remaining phrenic nerve, lung and diaphragm muscle function are sufficient to accommodate electrical stimulation. The candidate must be in otherwise generally good health.

The laparoscopic surgical procedure for implantation of the NeuRx Diaphragm Pacing System™ (DPS) is expected to take two hours. The patient is implanted with four IM (Intra-Muscular) electrodes, two in each hemidiaphragm, using laparoscopic techniques. Pre-operatives tests are obtained for patients based on their needs prior to general anesthesia.

14.2 ESTABLISH LAPAROSCOPIC PORTS

The operation is done in the supine position with no neuromuscular blocking agents. The patient's abdomen and chest is prepped and draped in the usual sterile fashion. Four ports will need to be inserted into the abdominal cavity: one for optics, two lateral working ports for the mapping probe and one epigastric port for electrode insertion instruments and the exit site from the abdominal cavity for the electrode lead wires. During this phase any abdominal adhesions are released and gastrostomy tube tracts are removed if they are in the way of implantation of the diaphragm pacing system. The falciform ligament is divided which allows easier visualization of the medial aspect of the right diaphragm and easier exit of the pacing electrodes through the epigastric port. Standard laparoscopic principles are followed with a typical setup.

4 ports insertion

15.0 MAPPING OF THE DIAPHRAGM

Prior to mapping, a surface electrode is placed around the thigh region and is connected to the Clinical Station (see Clinical Station Operating Manual). Intra-abdominal pressure during test stimulation (mapping) will be measured externally through the insufflation port with tubing. This tubing will connect to a pressure sensor to go to the clinical station. Mapping involves finding the point on the abdominal side of the diaphragm where stimulation causes the greatest diaphragm

excursion. The mapping instrument has a suction port that allows it to temporarily attach to the diaphragm and deliver an electrical stimulus. Stimulation is applied in either a twitch or burst mode from the clinical station. Mapping allows qualitative and quantitative data to be obtained. Quantitatively, changes in abdominal pressures are measured. Qualitatively, observation of the diaphragm contraction is performed. The stronger the stimulated contraction the closer to the motor point of the phrenic nerve. During mapping the entire diaphragm is assessed in grid pattern to be sure the point of maximal contraction is not missed.

These locations are recorded on transparency that is overlaid on the laparoscopic video display. The magnitude of the abdominal pressure change, to the applied stimuli, is recorded for each location. The primary electrode site is identified at the location of maximal pressure change in each hemi-diaphragm. A secondary electrode site is identified as either a backup to the primary site or at a location in each hemi-diaphragm that recruits another region (e.g. anterior, lateral, or posterior) of the diaphragm at a similar magnitude. The two locations are marked using a marker.

16.0 INSERTION AND PLACEMENT OF THE ELECTRODES

16.1 Loading Electrodes into the Electrode Delivery Instrument (Both instruments, Reusable or Disposable are loaded with the following technique)

Preparation of the Disposable Electrode Delivery Tools:

The scrub staff should load one Permaloc® Electrode into one of the Disposable Electrode Delivery Tools.

A bowl of sterile water or saline may be used to lubricate the electrodes prior to inserting into the Disposable Electrode Delivery Tool. Use a syringe filled with saline to lubricate the Disposable Electrode Delivery Tool to aid in electrode insertion. The Disposable Electrode Delivery Tool can be loaded with the needle locked at 90 or 180.


- Option #1: Loading electrode with needle at 180:
 1. Squeeze handle to the second locking position and lock the needle position by pushing the slide button.


2. Lay silver pin end of electrode into needle tip bevel area.
3. Keep electrode parallel to needle to not damage silicone sheath on electrode while loading through the tube of the tool.
4. DO NOT damage electrode by forcing through the needle bevel.
5. Continue to push electrode through the needle until the tip of the electrode nears the bevel of the needle.
6. Place thumb against bevel of the needle point and pull the remaining electrode through the tube of the tool. Ensure the blue tip is completely in the bevel of the needle.

7. Verify that only the electrode stimulating tip of the electrode is exiting the needle tip of the Disposable Electrode Delivery Tool.
8. While holding handle, unlock the needle position by pushing the slide button.
 - Option #2: Loading electrode with needle at 90:
 1. Squeeze handle to the first locking position and lock the needle by pushing the slide button.
 2. Lay silver pin end of electrode into needle tip bevel area.

3. Keep electrode parallel to needle to not damage silicone sheath on electrode.
4. DO NOT damage electrode by forcing through the needle bevel.
5. Push electrode through the needle until the tip of the electrode nears the bevel tip.
6. Insert silver pin end of the electrode into the electrode tube opening
7. Continue to push electrode through the electrode tube until exits tool.
8. Place thumb against bevel of the needle point and pull the remaining electrode through the tube of the tool. Ensure the blue tip is completely in the bevel of the needle.

9. Verify that only the electrode stimulating tip of the electrode is exiting the needle tip of the Disposable Electrode Delivery Tool.
10. Unlock the needle position by pushing the slide button.

Repeat this procedure for loading electrodes into other Disposable Electrode Delivery Tools.

16.2 Inserting Electrodes into the diaphragm using Electrode Delivery Instrument

Once the primary and secondary electrode sites are identified in each hemidiaphragm the implantation phase begins. An intramuscular electrode is introduced into the abdominal cavity with the electrode delivery instrument.

Surgical Electrode Delivery Instrument

The electrode is inserted into the diaphragm, at an angle so that the electrode lead travels parallel to the plane of the diaphragm prior to exit, and the delivery instrument is withdrawn. A standard laparoscopic dissector may be used to assist in the positioning of the electrode at the desired site on the diaphragm and give counter-traction when the needle is being removed. The electrode is then tested to assure the desired response to twitch stimuli is achieved and the procedure is repeated for the remaining electrodes. If the response is not adequate when tested the electrode may be withdrawn and another implanted. A second electrode is implanted at the previously marked site during mapping. At the conclusion of the implantation a CXR will be obtained to be sure no intra-abdominal air will have tracked with the needle to the chest cavity- a capnothorax. If one is present it can be aspirated with a percutaneous catheter (such as an angiocatheter, small chest tube, or thoracentesis) at the end of the case.


Electrode Insertion

Electrode Insertion

Once all four electrodes are implanted they are brought out through the epigastric trocar location.

17.0 ROUTING

The electrode leads are then routed to the percutaneous exit site. These lead wires will be tunneled subcutaneously to an area in the upper chest at a site deemed appropriate by the surgeon and patient's caregivers. An additional indifferent electrode will be placed subcutaneously in the upper chest through a separate percutaneous exit site. The electrodes are then retested to make sure that all of the connections have been made properly. An EKG strip is recorded with all four electrodes active to be sure there is no capture of the cardiac rhythm. At this time if the patient needs a gastrostomy a standard percutaneous endoscopic technique can be done and the wires are once again checked. The port incisions are closed and the patient is transferred to recovery. The wire exit site is covered with a wound dressing.

18.0 BLOCKING


Blocking is the final process of preparing the wires to allow connection to the NeuRx DPS and its connection cable. Blocking can be done in the Operating Room or in Recovery.

[CAUTION: When performing this step, avoid pulling excess electrode from the exit site. Excessive externalized electrode from the exit site could lead to accidental electrode breakage.]

1. If covered, carefully remove gauze pad to view the electrode exit site and clean the exit site with an alcohol pad. Obtain the crimping tool, the socket pusher and the electrode connector kit (22-0005).

Crimping Tool


Socket Pusher


Electrode Connector Kit


2. Remove a gold pin from the electrode connector kit. Using tweezers, grasp the larger end of the gold pin and place the tapered end into the crimping tool.


3. Using tweezers, gently insert the silver pin end of the electrode into the larger end of the gold pin that has been placed in the crimping tool.

4. While holding the silver pin into the gold pin, firmly squeeze the handle of the crimping tool until it stops. Release handle and carefully remove crimped electrode from crimping tool
5. Repeat this process until a gold pin has been crimped on all exiting electrodes.
6. Place the electrode connector into the socket pusher and orient the electrode connector with the two gold pins visible as pictured above.


7. Starting with the most cephalad (superior) electrode, insert the electrode into the 3rd hole in the electrode connector.

8. Gently close the socket pusher and line up the electrode into the slot of the pusher.
9. Once lined up, firmly close the socket pusher.



10. Select the next cephalad electrode and insert into the 5th hole in the electrode connector.



11. Gently close the socket pusher and line up the electrode into the slot of the pusher. Once lined up, firmly close the socket pusher.

12. Select the 3rd cephalad electrode and insert into the 6th hole in the electrode connector.


13. Gently close the socket pusher and line up the electrode into the slot of the pusher. Once lined up, firmly close the socket pusher.

14. Select the 4th cephalad electrode and insert into the 7th hole in the electrode connector.


15. Gently close the socket pusher and line up the electrode into the slot of the pusher. Once lined up, firmly close the socket pusher.

16. Select the indifferent electrode and insert into the 11th hole in the electrode connector.

17. Gently close the socket pusher and line up the electrode into the slot of the pusher. Once lined up, firmly close the socket pusher.


Assembled Electrode Connector


18. Carefully connect the patient cable to the electrode connector and perform a system check by using the NeuRx® Clinical Station at minimal settings. If an %+ is displayed in any channel, remove the identified gold pin from the electrode connector.

19.0 Apply Strain Relief Boot


1. The strain relief boot (found in the Surgical Connector Kit, part 22-0028) is pre-slit along the center of one of the flat surfaces to aid placement over the electrode wires.
2. Gently open the SLIT on the strain relief boot and wrap around the electrode leads and then slide the boot over the electrode connector block.

3. Turn the connector over (opposite of SLIT side) and hold the front flap of the boot up and remove the adhesive liner with forceps. Press the boot down firmly to the connector as shown.


4. Turn the boot back over to SLIT side. Remove the adhesive liner on the SLIT and firmly press the flap down. DO NOT DAMAGE THE ELECTRODES WHILE DOING THIS.


5. Remove the adhesive liner on the left side of the boot and press firmly to connector.


6. Remove the adhesive liner on the right side and press firmly on the connector as shown.


7. The assembled boot should appear as shown.



8. Remove the syringe filled with silicone adhesive from the kit and assemble grey tip to syringe. Save cap for resealing the tube after use.
9. Remove Plunger from bag and place in syringe. (Note: The Plunger is intentionally loose when placed in the syringe. It may fall out of syringe if turned upside down.)
10. Gently insert grey tip of syringe into strain-relief boot all the way to the electrode connector block.


11. While holding the strain relief boot, **gently** back-fill the strain relief boot from the bottom ensuring that there are no voids or air bubbles. Remove any excess silicone at the exit of the boot.

12. Remove backing from the Connector Holder and secure to patient.

13. Position the electrode connector block into the connector holder to dry and position the 5 wires into a flat and straight alignment.
14. Cover the exposed electrodes with a gauze pad and secure with a water-proof adhesive dressing.

20.0 REFERENCES

Onders RP, Aiyar H, Mortimer JT. Characterization of the Human Diaphragm Muscle with respect to the Phrenic Nerve Motor Points for Diaphragmatic Pacing. American Surgeon 2004;70:241-247.

Onders RP, Ignagni AI, Aiyer H, Mortimer JT. Mapping the Phrenic nerve Motor Point: the Key to a Successful Laparoscopic Diaphragm Pacing System in the first human series. Surgery 2004;136:819-826.

Onders RP, Ignagni AI, DeMarco AF, Mortimer JT. The learning Curve of investigational surgery: Lessons Learned from the first series of laparoscopic diaphragm pacing for chronic ventilator dependence. Surgical Endoscopy 2005;19: 633-7

A. J. Kerwin, B. K. Yorkgitis, D. J. Ebler, F. G. Madbak, A. T. Hsu, and M. L. Crandall, %Use of diaphragm pacing in the management of acute cervical spinal cord injury.,+J. Trauma Acute Care Surg., vol. 85, no. 5, pp. 928. 931, Nov. 2018.

R. P. Onders et al., %Completed FDA feasibility trial of surgically placed temporary diaphragm pacing electrodes: A promising option to prevent and treat respiratory failure.,+Am. J. Surg., vol. 215, no. 3, pp. 518. 521, Mar. 2018.

R. P. Onders, M. J. Elmo, C. Kaplan, R. Schilz, B. Katirji, and G. Tinkoff, %long-term experience with diaphragm pacing for traumatic spinal cord injury: Early implantation should be considered,+Surg. (United States), vol. 164, no. 4, pp. 705. 711, Oct. 2018.

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20.0 BIBLIOGRAFIA

Onders RP, Aiyar H, Mortimer JT. Characterization of the Human Diaphragm Muscle with respect to the Phrenic Nerve Motor Points for Diaphragmatic Pacing. American Surgeon 2004;70;241-247.

Onders RP, Ignagni AI, Aiyer H, Mortimer JT. Mapping the Phrenic nerve Motor Point: the Key to a Successful Laparoscopic Diaphragm Pacing System in the first human series. Surgery 2004;136;819-826.

Onders RP, Ignagni AI, DeMarco AF, Mortimer JT. The learning Curve of investigational surgery: Lessons Learned from the first series of laparoscopic diaphragm pacing for chronic ventilator dependence. Surgical Endoscopy 2005;19; 633-7

A. J. Kerwin, B. K. Yorkgitis, D. J. Ebler, F. G. Madbak, A. T. Hsu, and M. L. Crandall, %Use of diaphragm pacing in the management of acute cervical spinal cord injury.,+J. Trauma Acute Care Surg., vol. 85, no. 5, pp. 928. 931, Nov. 2018.

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R. P. Onders, M. J. Elmo, C. Kaplan, R. Schilz, B. Katirji, and G. Tinkoff, %Long-term experience with diaphragm pacing for traumatic spinal cord injury: Early implantation should be considered,+Surg. (United States), vol. 164, no. 4, pp. 705. 711, Oct. 2018.

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R. P. Onders et al., %Completed FDA feasibility trial of surgically placed temporary diaphragm pacing electrodes: A promising option to prevent and treat respiratory failure.,+Am. J. Surg., vol.215, no.3, pp.518. 521, Mar. 2018.

R.P.Onders, M. J.Elmo, C. Kaplan, R. Schilz, B. Katirji, and G. Tinkoff, %long-term experience with diaphragm pacing for traumatic spinal cord injury: Early implantation should be considered,+Surg. (United States), vol. 164, no.4, pp.705. 711, Oct. 2018.

 300 Artino Street
Oberlin, Ohio 44074
U.S.A.

www.synapsebiomedical.com

Tel: 1-888-767-3770
1-440-774-2488
Fax: 1-440-774-2572

Synapse Biomedical - Europe
7 Rue de la Liberation
95880 Enghien Les Bains, France

France and other Countries
Tel: +33 (0) 1 34 16 49 11
Fax: +33 (0) 1 74 18 08 19

© Synapse Biomedical Inc. 2021
All Rights Reserved
Tous droits réservés
Alle Rechte vorbehalten
Tutti i diritti riservati
Reservados todos los derechos

 77-0067 Rev G